



# BDO2 DRYING OVEN by BrugTech

## USAGE

### Efficient sample drying.

The BRUGTECH Drying Oven BDO2 is widely used in many laboratory and plant facilities. Samples are placed in trays and into one of the 2 steel trolleys.

Offering two recirculation fans for homogenous heat distribution and faster drying. Optionally, the drying oven can also be delivered with an added moisture damp extraction fan.

Fits up to  
144 samples

Up to 125° C  
Sample temp

## EASY OPERATION



## OPTIMIZE

- Provides the best drying conditions for all types of laboratory applications
- Uniform heat distribution
- Automatically control heating and timer functions
- Delivered with trolleys



BrugTech GmbH  
Rehagen 52  
30165 Hannover  
Germany



[info@brugtech.de](mailto:info@brugtech.de)



# BDO2 DRYING OVEN by BrugTech

## EASY OPERATION

## Safety

The **Oven** body is constructed with double hull steel plates filled in between with stone wool and glass wool for heat isolation purposes.

The inside of the **Oven** is covered with **stainless steel**. It is equipped with a movable trolley featuring 8 shelves - to easily load and unload samples. and avoid injuries.

**Safety Door Handle** - openable from the inside. Guaranteeing safety of your employees.

2-step heating limit for safety:

1. Digitally adjustable on the control panel
2. Analog heat limiter, avoids high temperature risks caused by operator. Oven will stop when the limit temperature is exceeded.



**BrugTech GmbH**  
Rehagen 52  
30165 Hannover  
Germany



[info@brugtech.de](mailto:info@brugtech.de)



# BDO2 DRYING OVEN by BrugTech



## EASY OPERATION



## Optimize

- Fits **2 x 8 shelves trolley**
- **Stainless Steel**
- Increased drying capacity with up to **9 trays per shelf**
- Allows **easy movement** of samples around the lab
- Easy to clean

BrugTech GmbH  
Rehagen 52  
30165 Hannover  
Germany



[info@brugtech.de](mailto:info@brugtech.de)